

**SP-1500-2 Calibration Chart**

\*Volume is the same in each compartment.

Height (in.)	Volume (Gal)						
61	1,685.82	50	1,381.82	38.75	1,070.91	27.5	760.00
60.75	1,678.91	49.75	1,374.91	38.5	1,064.00	27.25	753.09
60.5	1,672.00	49.5	1,368.00	38.25	1,057.09	27	746.18
60.25	1,665.09	49.25	1,361.09	38	1,050.18	26.75	739.27
60	1,658.18	49	1,354.18	37.75	1,043.27	26.5	732.36
59.75	1,651.27	48.75	1,347.27	37.5	1,036.36	26.25	725.45
59.5	1,644.36	48.5	1,340.36	37.25	1,029.45	26	718.55
59.25	1,637.45	48.25	1,333.45	37	1,022.55	25.75	711.64
59	1,630.55	48	1,326.55	36.75	1,015.64	25.5	704.73
58.75	1,623.64	47.75	1,319.64	36.5	1,008.73	25.25	697.82
58.5	1,616.73	47.5	1,312.73	36.25	1,001.82	25	690.91
58.25	1,609.82	47.25	1,305.82	36	994.91	24.75	684.00
58	1,602.91	47	1,298.91	35.75	988.00	24.5	677.09
57.75	1,596.00	46.75	1,292.00	35.5	981.09	24.25	670.18
57.5	1,589.09	46.5	1,285.09	35.25	974.18	24	663.27
57.25	1,582.18	46.25	1,278.18	35	967.27	23.75	656.36
57	1,575.27	46	1,271.27	34.75	960.36	23.5	649.45
56.75	1,568.36	45.75	1,264.36	34.5	953.45	23.25	642.55
56.5	1,561.45	45.5	1,257.45	34.25	946.55	23	635.64
56.25	1,554.55	45.25	1,250.55	34	939.64	22.75	628.73
56	1,547.64	45	1,243.64	33.75	932.73	22.5	621.82
55.75	1,540.73	44.75	1,236.73	33.5	925.82	22.25	614.91
55.5	1,533.82	44.5	1,229.82	33.25	918.91	22	608.00
55.25	1,526.91	44.25	1,222.91	33	912.00	21.75	601.09
55	1,520.00	44	1,216.00	32.75	905.09	21.5	594.18
54.75	1,513.09	43.75	1,209.09	32.5	898.18	21.25	587.27
54.5	1,506.18	43.5	1,202.18	32.25	891.27	21	580.36
54.25	1,499.27	43.25	1,195.27	32	884.36	20.75	573.45
54	1,492.36	43	1,188.36	31.75	877.45	20.5	566.55
53.75	1,485.45	42.75	1,181.45	31.5	870.55	20.25	559.64
53.5	1,478.55	42.5	1,174.55	31.25	863.64	20	552.73
53.25	1,471.64	42.25	1,167.64	31	856.73	19.75	545.82
53	1,464.73	42	1,160.73	30.75	849.82	19.5	538.91
52.75	1,457.82	41.75	1,153.82	30.5	842.91	19.25	532.00
52.5	1,450.91	41.5	1,146.91	30.25	836.00	19	525.09
52.25	1,444.00	41.25	1,140.00	30	829.09	18.75	518.18
52	1,437.09	41	1,133.09	29.75	822.18	18.5	511.27
51.75	1,430.18	40.75	1,126.18	29.5	815.27	18.25	504.36
51.5	1,423.27	40.5	1,119.27	29.25	808.36	18	497.45
51.25	1,416.36	40.25	1,112.36	29	801.45	17.75	490.55
51	1,409.45	40	1,105.45	28.75	794.55	17.5	483.64
50.75	1,402.55	39.75	1,098.55	28.5	787.64	17.25	476.73
50.5	1,395.64	39.5	1,091.64	28.25	780.73	17	469.82
50.25	1,388.73	39.25	1,084.73	28	773.82	16.75	462.91
		39	1,077.82	27.75	766.91	16.5	456.00

**SP-1500-2 Calibration Chart**

\*Volume is the same in each compartment.

16.25	449.09	5	138.18
16	442.18	4.75	131.27
15.75	435.27	4.5	124.36
15.5	428.36	4.25	117.45
15.25	421.45	4	110.55
15	414.55	3.75	103.64
14.75	407.64	3.5	96.73
14.5	400.73	3.25	89.82
14.25	393.82	3	82.91
14	386.91	2.75	76.00
13.75	380.00	2.5	69.09
13.5	373.09	2.25	62.18
13.25	366.18	2	55.27
13	359.27	1.75	48.36
12.75	352.36	1.5	41.45
12.5	345.45	1.25	34.55
12.25	338.55	1	27.64
12	331.64	0.75	20.73
11.75	324.73	0.5	13.82
11.5	317.82	0.25	6.91
11.25	310.91	0	0.00
11	304.00		
10.75	297.09		
10.5	290.18		
10.25	283.27		
10	276.36		
9.75	269.45		
9.5	262.55		
9.25	255.64		
9	248.73		
8.75	241.82		
8.5	234.91		
8.25	228.00		
8	221.09		
7.75	214.18		
7.5	207.27		
7.25	200.36		
7	193.45		
6.75	186.55		
6.5	179.64		
6.25	172.73		
6	165.82		
5.75	158.91		
5.5	152.00		
5.25	145.09		